



*Innovative, re-usable support for STAYONS
Leg Wrap application on knee & hock joints.*

Poultice Leg Wrap

Dry ingredient formula

- Lightweight, compact & easy to carry along
- Soak for just in 15-20 seconds, in warm / cold or ice water
- ***Dual Action***: Use with ice water to provide immediate cold therapy, then leave overnight for poultice benefits

Prefilled, preshaped

- Quilted pockets maintains poultice in close body contact
- Forms pliable poultice to hold moisture for extended time
- Molds beautifully around legs

Clean, rapid application & removal

- Poultice stays within wrap, no mess, saves time & effort

*Rapid, secure application & removal.
No mess. No problems.
4/pack - 12 packs/case*



Leg Wrap Application



1. SOAK. Using a STAYONS Poultrice Leg Wrap, fold wrap loosely and submerge in warm, cold or ice water as needed, for 15 to 20 seconds.

Poultrice rapidly absorbs and holds water. Remove wrap and allow excess water to drain back into pail.



2. WRAP. Wrap snugly around the leg, below knee.

STAYONS beautifully mold around the leg, providing maximum surface contact for the poultrice.



3. SECURE. Apply stable bandage to cover. *No brown paper or plastic wrap required. Saves extra work, materials, and time.*



DUAL-ACTION

In acute injury situations or to help reduce heat & inflammation from heavy training, soak STAYONS Leg Wraps in ice water. For maximum effect, soak wrap for 60 seconds.

The poultrice wrap will absorb & hold up to 3x its weight in ice water to provide immediate cold therapy benefits. As the heat from the leg dissipates, the clay/epsom poultrice blend continues to provide its drawing action, throughout the rest of the treatment period. Short term & long term relief for sore legs.

***Only STAYONS Leg Wraps provide this DUAL-ACTION support...
- and you only need to wrap the leg once!***