

# Telfa™ "Ouchless" Non-Adherent Dry Dressings



Order Info

Features

## Telfa™ "Ouchless" Non-Adherent Dry Dressings

Telfa™ "Ouchless" Adhesive Dressings

### Features and Benefits

Feature	Benefit
<i>Bonded on both sides with perforated non-adherent film</i>	Can be cut to fit any shape without separating
<i>Easy to use</i>	An ideal take home dressing
<i>Offers the most comprehensive range of sizes and offerings</i>	Can address most of the wound care and post surgical applications
<i>Superior "Ouchless" Telfa™ contact layer has clinically proven nonadherency</i>	Minimizes wound trauma; will not disrupt healing tissue by sticking to wound
<i>Polyester wound contact film polymer</i>	Can be steam sterilized
<i>Telfa™ film has 382 holes per square inch with a diameter of 0.020 – 0.030 inch</i>	Tested to provide the optimal balance of non-adherence and exudate migration
<i>Absorbent core is made out of 100% USP grade absorbent cotton</i>	Highly absorbent to help minimize skin maceration
<i>Ideal moisture vapor transmission rate</i>	Helps control the fluid loss and to provide moist wound healing environment

### Frequently Asked Questions

Question	Answer
<i>Why do Telfa™ Ouchless Non-adherent Dressings sometimes stick to the wound?</i>	Telfa™ dressings are designed for lightly exuding wounds. If the product is sticking, a more absorbent dressing may be indicated.
<i>What type of wounds are Telfa™ Ouchless Non-adherent Dressings designed for?</i>	This product is designed for light draining wounds such as surgical sites, skin tears or lightly draining wounds.