

SNAP* Bovine Virus Diarrhea Antigen Test Kit

IDEXX Labs

IDEXX SNAP* BVDV

Rapid in vitro diagnostic for the detection of bovine viral diarrhea virus in persistency-infected animals (BVD PI-animals) using ear-notch (ear tissue biopsy) and serum samples.

Note: for veterinary use only.

Kit Components:

One 42 mL bottle anti-BVDV: HRPO conjugate[†]

25 SNAP BVDV devices

Reagents contained in each device:

0.6 mL substrate solution

0.4 mL wash solution[‡]

Capped ear notch tubes

Transfer Pipettes

[†]Preserved with Bronidox and MIT

[‡]Preserved with ProClin™ 150

Storage:

2°C-8°C (35.6-46.4°F)