

J-LUBE OB POWDER

JORGENSEN LABORATORIES

A HAND LUBRICANT CONCENTRATE FOR VETERINARIAN USE ONLY

EASY TO USE

Apply to hands as a powdered hand soap. J-lube may be stirred into water to make a liquid lubricant for obstetrical use.

EASY TO CLEAN

J-lube is easily removed with water. A small amount of table salt added to the hands as a powdered hand soap will facilitate removal of all traces of the lubricant.

INGREDIENTS:

polyethylene polymer 25%
dispersing agent 75%

CAUTION: Leakage of this product into the peritoneal cavity may cause illness and/or death.