

NurseMate[®] ASAP with Immu-PRIME

Sterling Technology

Concentrated Bioactives for Newborn Calves

Use: To prime newborn calf's system to help assure maximum benefit from Mother's Colostrum.

Directions: Feed entire tube as soon after birth as possible by slowly placing it on the back of tongue. Feed at room temperature.

Ingredients: Bovine Colostrum Whey, Methyl Cellulose, Gum Blend, Sodium Benzoate (preservative), EDTA (preservative), Vitamin B-12, Vitamin H and Folic Acid.

Guaranteed Analysis: Each 30 ml syringe contains:

Crude Protein (Min.) .2%, Crude Fat (Min.) .1%, Crude Fiber (Max.) .1%

This product was made in a feed manufacturing facility that does not handle or store products containing animal proteins prohibited in ruminant feed.

Store in a cool, dry environment out of direct sunlight. Keep from freezing. Keep out of the reach of children.

Net Contents: 30ml