

Bovine IgG Colostrum 150 Replacer for Calves

OXFORD AG™

Made with whole bovine colostrum to feed newborn calves when maternal colostrum is unavailable or poor quality.

Colostrum 150 is made with whole bovine colostrum to provide your newborn calves with the immunity they need for a strong, healthy start!

Mixing Instructions

1. **FILL:** Remove cap from the Perfect Udder® Feeding System and add warm water (110-120 °F/ 43-49 °C) to the fill line. (1,875 ml)
2. **SHAKE:** Replace cap and shake and massage until all solids are dissolved completely.
3. **FEED:** Remove cap and attach either the nipple (Sold separately) or esophageal feed tube and be sure that the calf consumes the entire contents of the bag.

Indications

To replace maternal colostrum as an aid in the treatment of Failure to Passive Transfer (FPT) in newborn calves.

Contents

IgG ≥ 150g

IgG ≥ 26.6%

Feeding Instructions

Feed one bag within 2 hours of birth. Calves should receive a second feeding of colostrum within 8 hours.

Storage

Store out of Direct Sunlight in a dry location at or below 77 °F (25 °C).

Do not mix with other products.

A replacement dose is 1.2 bags.

For animal use only.

Product of Canada

Visit OxfordAg.com for video instructions on proper esophageal feed tube use.